

Introducing the Platinum Range

pelican
platinum



Comfort, security and flexibility.
So you can feel free to move in
your second skin.

pelican
HEALTHCARE



pelican
platinum

**The Platinum range
is all about comfort,
security and flexibility.**

It's your choice

"Pelican Platinum is far better than pouches I have tried before. I have confidence now, they are supple and stick better to my skin."

Platinum Contour Closed User

Closed range

VitaDerm Foam backed skin protector

- Advanced integrated hydrocolloid
- Adheres immediately to the skin and continues to improve
- Reduces leaks
- Increased security

Advanced filter system


- Controlled airflow filter
- Reduces risk of pancaking and ballooning
- Works as an excellent deodoriser
- Waterproof element reduces filter leakage

Contoured skin protector

- Features 6 flutes allowing for a more secure, snug fit
- Conforms to the natural body shape
- Ideal for surfaces with lumps, bumps or hernias

Soft and breathable fabric

- Non-woven fabric, which is comfortable and easy to dry.
- Durable yet breathable, minimising any embarrassing rustling sounds.
- Split fabric backing allows the pouch to be positioned more easily



Platinum Contour
Convex Closed

Platinum Closed

cut-to-fit
options
ranging from
12-80mm

Platinum Closed

	Size Guide	Clear	Opaque
Mini	Cut to fit 12-65mm	107212	107112
Standard	Cut to fit 12-65mm	107412	107312
	22.5mm	107422	107322
	25mm	107425	107325
	27.5mm	107427	107327
	30mm	107430	107330
	32.5mm	107432	107332
	35mm	107435	107335
	37.5mm	107437	107337
	40mm	107440	107340
	45mm	107445	107345
50mm	107450	107350	
55mm	107455	107355	
	Cut to fit 15-80mm	108415	108315
Maxi	Cut to fit 12-65mm	107612	107512

30 pouches per pack

Platinum Contour Closed

	Size Guide	Clear
Mini	Cut to fit 12-65mm	207212
Standard	Cut to fit 12-65mm	207412
	22.5mm	207422
	25mm	207425
	27.5mm	207427
	30mm	207430
	32.5mm	207432
	35mm	207435
	37.5mm	207437
	40mm	207440
	45mm	207445
50mm	207450	
55mm	207455	
	Cut to fit 15-80mm	208415
Maxi	Cut to fit 12-65mm	207612

30 pouches per pack

Platinum Convex Closed

	Size Guide	Clear	
Mini	Cut to fit 12-40mm	307212	
Standard	Cut to fit 12-40mm	307412	
	20mm	307420	
	25mm	307425	
	30mm	307430	
	35mm	307435	
	40mm	307440	
	Cut to fit 15-60mm	308415	
	Maxi	Cut to fit 12-40mm	307612

10 pouches per pack

Platinum Contour Convex Closed

	Size Guide	Clear	
Mini	Cut to fit 12-40mm	407212	
Standard	Cut to fit 12-40mm	407412	
	20mm	407420	
	25mm	407425	
	30mm	407430	
	35mm	407435	
	40mm	407440	
	Cut to fit 15-60mm	408415	
	Maxi	Cut to fit 12-40mm	407612

10 pouches per pack

Drainable range

VitaDerm Foam backed skin protector

- Advanced integrated hydrocolloid
- Adheres immediately to the skin and continues to improve
- Reduces leaks
- Increased security

Advanced filter system


- Controlled airflow filter
- Reduces risk of pancaking and ballooning
- Works as an excellent deodoriser
- Waterproof element reduces filter leakage

Contoured skin protector

- Features 6 flutes allowing for a more secure, snug fit
- Conforms to the natural body shape
- Ideal for surfaces with lumps, bumps or hernias

Discreet outlet

- Features a grey velcro tab, allowing the outlet to be folded away securely
- Centrally aligned outlet, suitable for use by colostomists and re-sited stomas



Platinum Contour
Convex Drainable



Platinum Drainable

cut-to-fit
options
ranging from
12-80mm

Platinum Drainable

	Size Guide	Clear
Mini	Cut to fit 12-65mm	118212
Standard	Cut to fit 12-65mm	118412
	20mm	118420
	22.5mm	118422
	25mm	118425
	27.5mm	118427
	30mm	118430
	32.5mm	118432
	35mm	118435
	37.5mm	118437
	40mm	118440
	45mm	118445
50mm	118440	
55mm	118455	
	Cut to fit 15-80mm	120415
Maxi	Cut to fit 12-65mm	118612

30 pouches per pack

Platinum Contour Drainable

	Size Guide	Clear
Mini	Cut to fit 12-65mm	218212
Standard	Cut to fit 12-65mm	218412
	20mm	218420
	22.5mm	218422
	25mm	218425
	27.5mm	218427
	30mm	218430
	32.5mm	218432
	35mm	218435
	37.5mm	218437
	40mm	218440
	45mm	218445
50mm	218450	
55mm	218455	
	Cut to fit 15-80mm	220415
Maxi	Cut to fit 12-65mm	218612

30 pouches per pack

Platinum Convex Drainable

	Size Guide	Clear
Mini	Cut to fit 12-40mm	318212
Standard	Cut to fit 12-40mm	318412
	20mm	318420
	25mm	318425
	30mm	318430
	35mm	318435
	40mm	318440
	Cut to fit 15-60mm	320415
Maxi	Cut to fit 12-40mm	318612

10 pouches per pack

Platinum Contour Convex Drainable

	Size Guide	Clear
Mini	Cut to fit 12-40mm	418212
Standard	Cut to fit 12-40mm	418412
	20mm	418420
	25mm	418425
	30mm	418430
	35mm	418435
	40mm	418440
	Cut to fit 15-60mm	420415
Maxi	Cut to fit 12-40mm	418612

10 pouches per pack

Urostomy range

VitaDerm Foam backed skin protector

- Advanced integrated hydrocolloid
- Adheres immediately to the skin and continues to improve
- Reduces leaks
- Increased security

Non-Return Valve

- Prevents the reflux of urine from travelling back up the pouch
- Reduces the risk of infection

Bung option

- Soft and flexible for a controlled flow
- Aids emptying with a universal section to a Night Drainage System
- Ideal for users with poor dexterity



Contoured skin protector


- Features 6 flutes allowing for a more secure, snug fit
- Conforms to the natural body shape
- Ideal for surfaces with lumps, bumps or hernias

Night Drainage Adaptor

- Soft and secure silicone adaptor
- Universal connection
- Hassle free – no need to change your regular night drainage bag
- Suitable for use with twist tap pouch

Twist-Tap option

- Secure tap features an "anti-drip" moulded stopper
- Red and green markers, allowing users to clearly see if the tap is **open** or **closed**
- Reduces the risk of accidents, removing any uncertainty for user



Platinum Contour
Convex Urostomy



Platinum Convex Urostomy

cut-to-fit
options
ranging from
12-60mm

	Size Guide	Platinum Convex Urostomy (Bung)	Platinum Contour Convex Urostomy (Bung)
Mini	Cut to fit 12-40mm	321212	421212
Standard	Cut to fit 12-40mm	321412	421412
	20mm	321420	421420
	25mm	321425	421425
	30mm	321430	421430
	35mm	321435	421435
	40mm	321440	421440
Maxi	Cut to fit 15-60mm	322415	422415
	Cut to fit 12-40mm	321612	421612

10 pouches per pack

	Size Guide	Platinum Convex Urostomy (Tap)	Platinum Contour Convex Urostomy (Tap)
Mini	Cut to fit 12-40mm	323212	423212
Standard	Cut to fit 12-40mm	323412	423412
	20mm	323420	423420
	25mm	323425	423425
	30mm	323430	423430
	35mm	323435	423435
	40mm	323440	423440
Maxi	Cut to fit 15-60mm	324415	424415
	Cut to fit 12-40mm	323612	423612

10 pouches per pack

The complete platinum range is now available in CONVEX


Platinum Convex

The Convex skin protector is soft, comfortable and flexible. This allows the user to forget they are even wearing a pouch and more importantly, reduces leakages by conforming to the body.

Suitable for flush or retracted stomas, ensuring a perfect fit.

All convex pouches should only be used after consulting a stoma care nurse.

pelican
platinum



"I found these pouches so comfortable and easy to use. They seemed to be more discreet than my current pouches."

Platinum Convex Drainable User



91%

of users rated Platinum as either EXCELLENT
or GOOD in terms of COMFORT.

It's your choice

*For complimentary samples
call us now on 0800 318 282*

 [pelicanhealthcareltd](https://www.facebook.com/pelicanhealthcareltd)  [@PelicanHealth](https://twitter.com/PelicanHealth)
[pelicanhealthcare.co.uk](https://www.pelicanhealthcare.co.uk)

pelican
HEALTHCARE