


an eakin company 

Paediatric & Neonatal

Designed with
children in mind

The Pelican Paediatric & Neonatal range offers an extensive range of pouches specifically designed for children.

0800 318 282

 pelicanhealthcareltd

pelicanhealthcare.co.uk

 @pelicanhealth

pelican
HEALTHCARE

Closed and Drainable



Skin Protector

The ultra-thin, hypoallergenic skin protector moulds to fit the body contours, providing a safe and secure fit.

Fabric

The soft fabric feels cool and comfortable against the skin, and is available in plain and teddy bear print.

Clipless System

Quick and easy to clean, the velcro tab under the fabric cover tucks the outlet away and holds it securely in place.



Waterproof Filter

The Drainable version features an improved filter system with a controlled flow for greater comfort. This feature significantly reduces the risk of pouch ballooning and pancaking.

Urostomy



Shape

The advanced pouch shape which is designed to minimise pooling, incorporates a comfortable backing behind the bung to ensure comfort.

Bung

The outlet features a soft and flexible bung, and by squeezing, gives you the ability to control flow when emptying.

*For complimentary samples,
call us now on 0800 318 282*

pelican
HEALTHCARE